

# ARTISTRY LABS

## AQUABRASION DEVICE

Multi-function skincare device that maximizes at-home skincare experience with 2 separate modes.

**97%** of clinical study panelists showed improvement in appearance of fine lines.\*

**86%** of consumers agreed that they prefer doing ARTISTRY LABS System treatment than visit a dermatologist.\*

\*based on clinical study and consumer test




### Suction Mode

Exfoliates loosened dead skin cells and extracts pores for deep cleansing  
Comes with 2 different sized heads for ease of fit on different parts of your face

### Serum Mode


Boosts the activity of key ingredients in Retexturizing Serum up to 900%\* or 10x\* more antioxidant and skin barrier functions

\*in vitro study with lentil seed extract.



# 7-DAY SKINCARE ROUTINE

Repeat this weekly routine for one month



Cleanser

Toner

Retexturizing Serum

Moisturiser

SPF



DAY	1	2 3	4	5 6 7
PEEL DAY				



## SUCTION Mode Treatment Instructions:

Avoid the eye area

\*Use your free hand to hold skin firm

Large Nozzle

- Forehead
- Cheeks
- Chin area
- Mouth area
- Nose area

Small Nozzle

## SERUM Mode Treatment Instructions:

Avoid the eye area

Avoid the centre of the neck area

Lightly touch the Dome on the face

Serum Dome

3 seconds beep

1 inch

3 seconds beep

1 inch

3 seconds beep

- Forehead
- Cheeks
- Chin area
- Nose area
- Neck area